


Reference exhibition

3D Gel Dosimeter for Radiotherapy

A polymer gel dosimeter is one of the chemical dosimeters and can be read out with high accuracy using MRI.




X-ray irradiation



3D Dose distribution

- **True 3D dosimetry**
 - No directional dependence
 - Homogeneous gel matrix containing reactive monomers
- **Room temperature storage**
 - This is the first polymer gel dosimeter that can be stored at room temperature.
 - A thermally tolerant gel matrix and a chemically stable oxygen scavenger are used.
- **Tissue equivalence**
 - Main component is water.

● Comparison of dose distribution




● Dosimetry system

- A pair of gel dosimeter bottles: one for preparing the dose-response curve, and another for dose verification
- MRI image processing software
- Solid phantom (optional)


● Chemical compositions of the gel dosimeter

- Reactive monomers
- Hydrogel: no gel-to-sol transition upon heating
- Oxygen scavenger: stable at room temperature


- **Dose-response curve**



- **No dose-rate dependence**



- **No integral dose dependence**



● **Measurement steps**

1. **X-ray irradiation**: Generation of hydroxyl radical resulting from water radiolysis
2. **Chemical reaction**: Polymerizing a reactive monomer and a cross-linker in the hydrogel induced by the hydroxyl radical
3. **MRI measurement**: Measuring the T2 relaxation time with MRI
4. **Image processing**: MRI image processing with the software converting the T2 to the corresponding relaxation rate R2

● **References of polymer gel dosimetry**

Polymer gel dosimetry

L J Schreiner, T Olding and K B McAuley

2010 *J. Phys.: Conf. Ser.* **250** 012014

Topical Review: Polymer gel dosimetry

C Baldock, Y De Deene, S Doran, G Ibbott, A Jirasek, M Lepage, K B McAuley, M Oldham, and L J Schreiner

2010 *Phys Med Biol.* **Feb 11**; 55(5): R1–R63

Polymer gel – TPS radiotherapy dosimetry GeVero® software for ionizing radiation absorbed dose 3D distribution calculations

Marek Kozicki, Piotr Maras, Jacek Jankowski and Andrzej C Karwowski

2009 *J. Phys.: Conf. Ser.* **164** 012062

Data from RIKEN

July 12, 2023

 **KYOTO KAGAKU co.,LTD**

Worldwide Inquiries & Ordering

<http://www.kyotokagaku.com>

E-mail: rw-kyoto@kyotokagaku.co.jp

Kyotokagaku Head Office and Factories:

15 Kitanechoya-cho, Fushimi-ku, Kyoto, 612-8388, JAPAN

Tel: +81-75-605-2510 Fax: +81-75-605-2519

Kyotokagaku USA Office:

USA, Canada and Mexico sales and services

3109 Lomita Boulevard, Torrance, CA 90505-5108, USA

Tel: 1-310-325-8860 Fax: 1-310-325-8867