



# Cardiology Patient Simulator K Version II


## Hone your senses and skills for bedside cardiology

See, listen, and feel

The cardiology patient simulator K II provides comprehensive hands-on training in cardiology and its complications

### Product Supervision:

Japanese Educational Clinical Cardiology Society,  
Toshio Amano, Masaya Kino, Takaharu Saito, Tsunekazu Takashina



### CASES

(MODE 1) Comprehensive patient simulation with sounds, pulses, apex beats and ECG

| No.  | Normal heart simulation (12 cases) | No.  | Heart disease simulation (14 cases) | No.  | Arrhythmia (10 cases)                 |
|------|------------------------------------|------|-------------------------------------|------|---------------------------------------|
| A-01 | S2 split (-) HR: 60                | B-01 | aortic stenosis                     | C-01 | sinus arrhythmia                      |
| A-02 | S1 split (+)                       | B-02 | mitral regurgitation                | C-02 | sinus tachycardia                     |
| A-03 | S2 split (+)                       | B-03 | mitral stenosis                     | C-03 | sinus bradycardia                     |
| A-04 | S2 wide split                      | B-04 | aortic regurgitation                | C-04 | ventricular premature contraction (1) |
| A-05 | S3 gallop                          | B-05 | hypertrophic cardiomyopathy         | C-05 | ventricular premature contraction (2) |
| A-06 | S4 gallop                          | B-06 | mitral steno-regurgitation          | C-06 | ventricular premature contraction (3) |
| A-07 | pulmonic ejection sound            | B-07 | pulmonic valvular stenosis          | C-07 | sino-atrial block                     |
| A-08 | S3 and S4 gallop                   | B-08 | atrial septal defect                | C-08 | atrio-ventricular block               |
| A-09 | innocent murmur                    | B-09 | ventricular septal defect           | C-09 | atrial fibrillation                   |
| A-10 | midsystolic click sound            | B-10 | tricuspid regurgitation             | C-10 | atrial flutter                        |
| A-11 | S2 split (-) HR: 72                | B-11 | acute mitral regurgitation          |      |                                       |
| A-12 | S2 split (-) HR: 84                | B-12 | patent ductus arteriosus            |      |                                       |
|      |                                    | B-13 | mitral valvular prolapse            |      |                                       |
|      |                                    | B-14 | dilated cardiomyopathy              |      |                                       |


(MODE 2) Arrhythmia simulation, auscultation training with ECG

| A                        | B                        | C                     | D                        |
|--------------------------|--------------------------|-----------------------|--------------------------|
| A-01 normal sinus R      | B-01 atrial flutter      | C-01 VVI pacemaker    | D-01 vpc (quadrigeminy)  |
| A-02 sinus tachycardia   | B-02 av block            | C-02 atrial pacemaker | D-02 vpc (trigeminy)     |
| A-03 sinus arrhythmia    | B-03 av block & crbbb    | C-03 vent pacemaker   | D-03 vpc (bigeminy)      |
| A-04 apc solitary        | B-04 av block (digial)   | C-04 av seq pacemaker | D-04 vpc (couplet)       |
| A-05 apc bigeminy        | B-05 av block (mobitz)   | C-05 icrbbb           | D-05 pvc (repetitive)    |
| A-06 ectopic pacemaker   | B-06 av block (mobitz)   | C-06 crbbb            | D-06 pvc (R-on-T type)   |
| A-07 wondering pacemaker | B-07 av block (3:1&4:1)  | C-07 clbbb            | D-07 non-sustained VT    |
| A-08 coronary sinus R    | B-08 av & crbbb          | C-08 clbbb            | D-08 vent tachycardia    |
| A-09 sinus bradycardia   | B-09 paroxysm atr tachy  | C-09 clbbb (by ami)   | D-09 vent flutter        |
| A-10 ss syndrome         | B-10 av junc R(svt)      | C-10 wpw syndrome     | D-10 vent fibrillation   |
| A-11 atrial fibrillation | B-11 av junc R (pat)     | C-11 wpw syndrome     | D-11 vent R (sinus cond) |
| A-12 atrial flutter      | B-12 av junc R           | C-12 wpw syndrome     | D-12 accel vent rhythm   |
| A-13 atrial flutter fib  | B-13 ac junc contraction | C-13 vpc (solitary)   | D-13 agonal rhythm       |

## Basic Skill Training

## Physical findings of Simulator K Version II

Comprehensive Clinical Examination Training



## Heart Sounds and Murmurs

Auscultation synchronized with electrocardiograms and findings

- differentiation of S1, S2, S3 and S4
- distinction between systolic and diastolic murmurs

## Pulse examination

Slight variations of the arterial pulse waves under different cardiac conditions or arrhythmias can be detected by pulse examination

## Observation of Jugular Veins

The strength and timing of "α" and "V" waves vary with cases

## Palpation of Apex Beat

Simulated RV, LV and DLV can be palpated

## Monitoring Screen

Real time dynamic charts on the screen:

- phonocardiogram
- sphygmogram (CAP/JVP/AC)
- electrocardiogram (ECG)

## Respiratory Sounds and Abdominal Movement

Respiratory cycles are represented to understand respiratory conditions

## For Effective Training Sessions

## Up to 5 Simulators at One Time

OSCE, group session, and scenario-based simulation



Up to five cardiology simulated cases can be controlled by one wireless tablet. Different conditions can be set in each simulator

The system allows for:

- scenario-based training involving more than one patient
- skills assessment sessions with different cardiovascular stations

## Playlist System

Simulator "K" is not just for individual task training

The incorporated "Play List Maker" facilitates creation of scenario-based training sessions with conditions that vary over time

## Playlist Maker facilitates:

- sessions with temporal change in physical findings
- standardization of training contents among different instructors
- saving time for preparation

## Ready-to-use

Simple connection and all-in-one unit structure

## Error Indicator

The error indicator performs check-up of the system to keep Simulator K in its best condition

## DESCRIPTIONS

## SPECIFICATIONS

Size: W145 x D75 x H170 cm / W57 x D30 x H70 in.

Weight: 145kg/320lbs

Power: AC100V-240V 50/60Hz

Power Consumption: 400W

## MATERIALS

Soft resin, hard resin  
Latex free

## SET INCLUDES

- 1 main unit with the model body
- 1 external speaker (built-in)
- 1 PC / keyboard / mouse (built-in)
- 1 air compressor unit (built-in)
- 1 main monitor

- 1 tablet
- 1 rib sheet
- 1 storage cover
- 4 textbooks
- 2 instruction manuals
  - built-in help file of software operation
  - instruction manual

\* A stethoscope is not included.