



Pitting Edema Model

Five variations of pitting edema models including a normal one. Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"

The models can be worn on the lower thighs of a SP or a training manikin.

Training items:

- Differentiation of pitting edemas with four different grades and a normal.
- Learn the examination procedures for edemas on the lower thighs.



Production Supervision:

Junji Otaki

ProfessorTokyo Medical University

Features:

- Four grades of pitting edemas with different resistances, depth of indentations.
- Grades of the edemas are based on “Mosby,s Guide to Physical Examination (7th Ed.)”
- The models can be worn on the lower thighs of a SP or a training manikin.
- The removable fibula plate facilitates training without SP leg or simulation on the soft bases.



Degrees of pitting edemas:



Normal
No indentation left



+ Light pitting
Indentations with approx. 2mm depth



2+ Moderate pitting
Indentations with approx. 4mm depth



3+ Deep pitting
Veins, bones and joints are indistinct because of swell.
Indentations with approx. 6mm depth.




4+ Very deep pitting
Indentations with approx. 8mm depth.

Set includes:

Strap-on pitting edema model 5 variations — 1 piece each



 **KYOTO KAGAKU** www.kyotokagaku.com rw-kyoto@kyotokagaku.co.jp

Head Office and Factories:
15 Kitanechoya-cho, Fushimi-ku, Kyoto, 612-8388, JAPAN
Tel: +81-75-605-2510 Fax: +81-75-605-2519

USA Office:
3109 Lomita Boulevard, Torrance, CA 90505-5108, USA
Tel: 1-310-325-8860 Fax: 1-310-325-8867