


41350-200-16 (BMI 32)/ 41350-200-17 (BMI 40)

Optional Parts for PH-2/ 2B/ 2E/ 60

Body plates

Body plates to simulate a patient of BMI 32/ BMI 40



FEATURES

For study of effect of patient size on radiation dose and image quality.

DESCRIPTIONS

SET INCLUDES

- 1 body plate (front)
- 1 body plate (back)
- 2 belts

MATERIALS

Urethane based resin (specific gravity: 1.06)

SPECIFICATIONS

Phantom size:
W41 x D45 x H28cm
W16.1 x D17.7 x H11 in

Packing size:
77 x 48 x 40 cm x 2 boxes
30.3 x 18.9 x 15.7 in x 2 boxes


Packing weight:
41350-200-16(BMI 32): 21 kg/ 46.3 lb
41350-200-17(BMI 40): 34.5 kg/ 76 lb

41350-000-11

Optional Parts for PH-2/ 2B/ 2E/ 60

Fractured Hand/Forearm Phantom PH-2/ 2B

X-ray phantom for trauma evaluation



DESCRIPTIONS

Bone Fractures:

ulna, radius, first metacarpal, middle phalanx of the index finger, distal phalanx of the first finger (compressed fracture), fifth metacarpal

SET INCLUDES

- 1 fractured hand/forearm phantom

MATERIALS

Soft tissue: urethane based resin (specific gravity: 1.06)
Synthetic bone: epoxy resin (specific gravity: 1.31)
*Phantom has no metal parts or liquid structure

41363-070

Optional Parts for PH-2/ 2B/ 2E/ 60

Storage case 2 (a pair)

