




Acquire skills for safe pediatric LP

Simulator for acquiring safe and reliable skills for pediatric lumbar puncture

Features

1. Soft, seamless baby manikin facilitates trainings in safe patient positioning and holding.
2. Two types of puncture pads with deep and shallow depth.
3. Realistic resistance to needle tips
4. Anatomically correct landmarks provides excellent hands-on experience.

Level of posterior superior iliac crests



Product Supervision

Takahiro Amano, M.D.
Vice President
Senior Vice-Dean, Postgraduate School
Professor and Director, Center of Postgraduate
Medical Education International University of Health and Welfare
Honorary Director, Sanno Medical Center

Gregory A. Plotnikoff, M.D., M.T.S., F.A.C.P.
Medical Director, Abbott Northwestern Hospital

Training Skills

- Patient positioning
- Palpation of landmarks
- Lumbar puncture and CSF collection
- CSF pressure measurement



Palpation of landmarks



Lumbar puncture and CSF collection



CSF pressure measurement

Specifications :

Set Includes:		Manikin Size:	Replacement Parts:
1 baby manikin (7-10 months)	1 stand for irrigation bag	W43 x D34 x H20 cm	(11348-410)
2 pediatric LP blocks (shallow)	1 syringe (50mL)		2 pediatric LP blocks (shallow)
2 pediatric LP blocks (deep)	1 instruction manual		spinal cord at 12 mm depth
1 irrigation bag			(11348-420)
			2 pediatric LP blocks(deep)
			spinal cord at 17 mm depth

* Needles are not included.