

M50B-B



Simulator Intravenous Arm II



Set of two models for variable practice

Palpable veins to simulate venipuncture, intravenous cannulation and injection.

FEATURES

- 1 | 2 kinds of puncture pads for difficulty variation
- 2 | 3 accessible venous lines at cubital fossa: basilic vein, median vein, and cephalic vein
- 3 | Realistic needle tip resistance of the vein walls
- 4 | "Flash-back" can be observed.

SKILLS

- | Intravenous injection, blood collection and drip infusion

DESCRIPTIONS

SPECIFICATIONS

Size: W54 x D11 x H11 cm each
W21.3 x D4.3 x H4.3 inch

Weight: kg/ lbs

Power: AC100V, 50/60 Hz
Power Consumption:

SET INCLUDES

- | | |
|-------------------------------------|----------------------------------|
| 2 arm models with supporter stand | 1 set of spoon and plastic jar |
| 2 injection pads | 1 instruction manual |
| 2 injection pads (red) | |
| 1 circulation pump | |
| 1 jar of coloring powder (red) | |

MATERIALS

Soft resin
Latex free

RECOMMENDED DEVICES

-

REPLACEMENT PARTS

- | | |
|--------------|--|
| 11267-300-03 | 2 replacement injection pad for M50B |
| 11267-300-04 | 2 replacement injection pad for M50B (red) |
| 11266-080 | 1 jar of coloring powder |

M50D

Basic IV Arm

Quality and economy in one

Basic IV Arm is realistic with fine detail, inexpensive and durable, perfect for a large number of trainees. Lightweight, portable and easy set up so that each trainee can have their own arm to practice. The vein tube can easily be filled with simulated blood using a soft plastic injection bottle.



FEATURES

- 1 | Flexible skin and palpable veins (replicates venous dilatation)
- 2 | Realistic needle tip sensation and flashback
- 3 | Injection pad is easily replaceable and durable.
- 4 | Simulate hemostasis by compressing the injection site
- 5 | Adhesive plaster sticks well to the arm and can be easily removed.
- 6 | Integrated stand for ease of set up
- 7 | Compact, lightweight, and portable

SKILLS

- | Tourniquet application
- | Venipuncture (including puncture with IV cannula)

DESCRIPTIONS

SPECIFICATIONS

Size: W54 x D11 x D11 cm
W21.3 x D4.3 x H4.3 inch

Weight: kg/ lbs

Power: -
Power Consumption: -

SET INCLUDES

- | | |
|------------------------------------|--|
| 1 arm model with supporter stand | 1 injection bottle for simulated blood |
| 1 injection pad | 1 carrying bag |
| 1 syringe (50ml) | |
| 1 jar of coloring powder (red) | |
| 1 spoon | |

MATERIALS

Soft resin
Latex free

RECOMMENDED DEVICES

-

REPLACEMENT PARTS

- | | |
|-----------|--------------------------|
| 11388-200 | 2 injection pads |
| 11266-080 | 1 jar of coloring powder |

