



Five variations of pitting edema

Features

1. Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"
2. The models can be attached to the lower thighs of an SP or a training manikin.
3. Pitting edema is classified into 4 levels by the degree of response to pressure and the depth of indentations.
4. The detachable fibula plate facilitates training simulation on a softer base.

Training Skills

- Differentiation of the pitting edema: four different levels and a normal case
- Examination procedures for edemas on the lower thighs

Product Supervision

Junji Otaki
Professor of Tokyo Medical University



Normal

1+ Light Pitting

2+ Moderate Pitting



3+ Deep Pitting

4+ Very Deep Pitting

Specifications :

Set Includes:	Manikin Size:
5 variations of strap-on pitting edema model 1 piece each	W10 x D23 cm each



Four stages of pressure sores for training

Features

1. The set consists of four types of models that explain each case and four types of sores to put on a patient care manikin.
2. The set facilitates both the classroom lecture and the practice according to each stage of the pressure sore.
3. Water color is included in the set.

Training Skills

- Pressure sores care



Specifications :

Set Includes:	Manikin Size:
4 types of the models for explanation	W38 x D38 cm each
4 types of the sore skins model	
1 set of water color paints	
1 set of brushes	
1 set of palettes	