


MW20

# Prostate and Rectal Simulator

## 2-in-1 simulator for prostate and rectal exams

The Prostate and Rectal Simulator includes 9 types of prostates and 4 types of rectums which can easily be changed. The lower torso manikin can be positioned in 3 different ways



### FEATURES

- 1 | Nine types of prostate simulate different scenarios
- 2 | Four rectum units
- 3 | Training with multiple prostates and rectal cases are concurrent and time efficient
- 4 | Lateral, supine and prone positioning
- 5 | 3 revolving units with 9 prostates can easily be changed.

### SKILLS

- | Digital-examination of the prostate and rectum
- | Insertion and the use of the anal-scope and the proctoscope


### KEY FEATURES

#### Cases: Nine types of Prostates



- |                  |                |                |
|------------------|----------------|----------------|
| #1 normal        | #4 prostatitis | #7 carcinoma 3 |
| #2 enlargement 1 | #5 carcinoma 1 | #8 carcinoma 4 |
| #3 enlargement 2 | #6 carcinoma 2 | #9 carcinoma 5 |

#### Cases: Four Rectum Units



- |                    |  |
|--------------------|--|
| #1 normal          | #3 large carcinoma                     |
| #2 small carcinoma | #4 polyp combines with small carcinoma |

### DESCRIPTIONS

#### SPECIFICATIONS

Size: 34x33x24cm/13.4x13x9.4 inch

#### SET INCLUDES

1	lower torso manikin	4	rectum units
3	prostate units	1	Vaseline
	(include 9 types of prostates)	1	body
		1	storage case

#### MATERIALS

Soft resin  
Latex free

#### REPLACEMENT PARTS

11399-010 prostate unit A(normal, enlargement 1&2)  
11399-020 prostate unit B(prostatitis, carcinoma 1&2)  
11399-030 prostate unit C(carcinoma 3,4&5)  
11399-040 rectum (normal)  
11399-050 rectum (small carcinoma)  
11399-060 rectum (large carcinoma)  
11399-070 rectum (polyp and small carcinoma)  
11399-080 skin