

MW22

# Patient Care Simulator "Cherry"

## Basic patient care manikin for introductory nursing

A basic full-body, nursing manikin that offers trainings of patient care for health care professionals



### FEATURES

- 1 | Able to sit in upright position
- 2 | Seamless limbs up to fingers and toes
- 3 | TPN and IV routes to simulate patients with catheter lines
- 4 | True-to-life articulation for patient handling




Subcutaneous injection



Wound care

### SKILLS



### DESCRIPTIONS

#### SPECIFICATIONS

Size: 157 cm/62 inch

#### SET INCLUDES

1	manikin
1	denture
1	injection pads (subcutaneous/ gluteal/ femoral IM/ Upper arm IM)

#### MATERIALS

Soft resin  
Latex free

1	Female genitalia
1	Male genitalia
1	Wound care pad
1	Bottle
1	Syringe

#### REPLACEMENT PARTS

11401-010	Insulin injection pad for MW22 (a set of 6)
11401-020	Gluteal IM pads for MW22 (a set of 4)
11401-030	Female genitalia unit for MW22
11401-040	Male genitalia unit for MW22
11401-050	Shoulder IM pads for MW 22 (a pair)
11401-060	Thigh IM pads for MW22
11401-070	Wound care pad for MW22 (a pair)
11401-080	Dentures(a set of upper & lower) for MW22, MW26
11401-090	Stomach for MW22
11388-200	Median antebrachial pad (a pair)