

MW2A



Male Catheterization/Enema Simulator

MW2B



Female Catheterization/Enema Simulator



MW2A



MW2B

High quality simulators for urethral catheterization, perineal care and enema practice

Achieve better training results with our soft and realistic catheterization simulator that can also be used on SP or full body manikins.

FEATURES

- 1 | Soft and flexible genitalia allow proper cleaning and procedural preparation
- 2 | Successful insertion confirmation by fluid outflow
- 3 | Enema insertion in lateral position
- 4 | Detachable external genitalia

SKILLS

- | | |
|---------------------------------|----------------------------------|
| Intermittent catheterization | Cleansing and aseptic technique |
| Indwelling catheterization | Catheter insertion and placement |
| Perineal care | Catheter Removal |
| Enema | Manual bladder compression |
| Localization of urethral meatus | |

Compact sets consist of a genitalia unit and accessories, for hybrid training and table top demonstration.

MW2C 11381-199-B

Male Genital Unit only

**MW2D** 11381-299-B

Female Genital Unit only



Enema Training



Pass of catheter can be observed through transparent bladder reservoir

**DESCRIPTIONS****SPECIFICATIONS**

Size: 34x33x24cm/ 13.5x13x9.5 inch

SET INCLUDES

- | |
|--|
| 1 lower torso manikin |
| 1 genital organ unit |
| 3 valves (one is incorporated in the assembled unit) |
| 1 irrigation bottle |

MATERIALS

Soft resin
Latex free

- | |
|-----------------------------|
| 1 rectum tube |
| 3 irrigation drainage tubes |
| 1 lubricant |
| 1 storage case |
| 1 instruction manual |

REPLACEMENT PARTS

- | | |
|--------------|---|
| 11381-100-01 | clip (small) for urethral tube for MW2A |
| 11381-100-02 | clip (large) for rectum unit |
| 11381-100-03 | urinary bladder reservoir |
| 11381-100-04 | rectum unit |
| 11381-300 | 6 replacement valves for MW2A/B |
| 11381-610 | 1 replacement female skin |
| 11381-700 | 2 replacement female meatus |
| 11381-800 | lubricant |
| 11381-410 | 1 replacement male skin |
| 11381-510 | 2 replacement male meatus |

RECOMMENDED DEVICES

* Catheters are not included.