MW9

Peripheral Venous Catheter Placement Simulator

Holistic simulation of peripheral intravenous access management in just one arm



FEATURES

- 1 | IVC insertion and management from start to finish
- 2 | Two puncture sites: median antebrachial vein and dorsal veins
- 3 | Durable puncture pad

SKILLS

- Infection control
- | Use of tourniquet
- | Identifying puncture site
- Sterilization of puncture site
- Puncture with IV cannula
- | Confirmation of flashback
- Application of pressure on vein
- Preparing and administering IV medications (simulated)
- | Dressing and securing cannula



Confirmation of flash back



Fixation



Application of pressure

DESCRIPTIONS

SPECIFICATIONS

SET INCLUDES

Size: W52 x D10 x H10 cm W20.5 x D3.9 x H3.9 inch

Weight: kg/ lbs

Power: AC100V-240V, 50/60 Hz Soft resin Power Consumption:

MATERIALS Latex free

10 simulated blood sticks (swab type)

- carrying bag
- 1 instruction manual


RECOMMENDED DEVICES

REPLACEMENT PARTS


11388-200 2 injection pads

for median antebrachial vein 11388-300 2 injection pads

for dorsal veins of hand 11388-400 10 simulated blood sticks (swab type)



1 bottle for simulated blood



infusion set

infusion bag

syringe (50ml) 1 plastic beaker



