





**Durable and stable**

Useful both for daily assessment and further research. Gray scale for evaluation, cyst targets with non-resonance cylinders, line targets for geometrical evaluation including close range (dead zone) resolution, axial and lateral resolution are prepared for scanning. The phantom is designed to allow scanning from all four side walls.



**Specifications :**


Set Includes:	Phantom Size:
1 phantom	W29 x D7 x H22 cm
1 storage case	





**Ensure highly detailed images to ensure reliable breast cancer examinations**

Ultrasound QA phantom for high precision imaging in the high frequency sonography around 10 MHz that is required in breast examination. For monthly basic quality check of ultrasound images, as well as longer term quality assurance to maintain consistency of the performance of scanners and transducers.

**Mass targets block (contrast resolution)**



**Mass targets block (contrast resolution)**



**Specifications :**

Set Includes:	Phantom Size:
1 mass targets block	W18 x D7.5 x H11 cm
1 dot targets block	
1 thermometer	
1 storage case	

\* Ultrasound machines is not included.