


MW50

Strap-on Venipuncture Simulator

Cost-effective training unit for venipuncture

Strap-on Venipuncture Simulator is a simple, economical training kit for intravenous injection. This trainer is to be strapped onto the upper arm of an SP (Standardized Patient) or a life-size training manikin.



FEATURES

- 1 | Realistic tissue and vein wall resistance
- 2 | Vein tube can be filled and refilled easily using a soft plastic bottle that works as a manual pump.
- 3 | Two vein tubes are embedded.

SKILLS

- | Blood collection and intravenous injection

DESCRIPTIONS

SPECIFICATIONS

Size: W12 x D14 x H2.5 cm / W4.7 x D5.5 x H0.98 inch
W6.5 x D11 x H2.0 cm / W2.6 x D4.3 x H0.79 inch (Pad)

SET INCLUDES

- 5 puncture pads (vein tube embedded)
- 5 supporting frames
- 5 belts
- 5 injection bottles
- 1 instruction manual
- 1 storage case
- 5 syringes

MATERIALS

Soft resin
Latex free

REPLACEMENT PARTS

11431-010 5 replacement puncture pads

OPTIONAL PARTS

11430-010 Arm Vessels Anatomy Sleevelet
11388-400 10 cotton swab red

RECOMMENDED DEVICES



-

M94

Intradermal Injection Training Simulator

Wheal can be shown by correct injection

This simulator provides training for intradermal injection to the inner forearm. Swelling of the skin will form when proper injection procedures are followed.



FEATURES

- 1 | Life like touch and resistance of the injection site resembles the human arm
- 2 | Injection pads are easily changeable.
- 3 | No needle trace remains.

SKILLS

- | Intradermal Injection

DESCRIPTIONS

SPECIFICATIONS

Size: W9 x D20 x H8 cm Weight: 0.3 kg/ 0.7 lbs
W3.5 x D7.9 x H3.1 inch

SET INCLUDES

- 5 injection pads
- 1 forearm model
- 1 supporter stand
- 4 injection pad cases
- 1 instruction manual

MATERIALS

Soft resin
Latex free

REPLACEMENT PARTS

11346-010 1 replacement puncture pad

RECOMMENDED DEVICES

Needle 29G