

PH-17 | 41317-000

## Water Body Phantom WAC



### Image evaluation in plain X-ray



SHOW MORE!



#### FEATURES

Water body phantom represents human chest and abdomen to serve as radiation absorber and scatterer.

#### DESCRIPTIONS

##### SET INCLUDES

- 1 body phantom
- 1 storage case

##### SPECIFICATIONS


Phantom size:  
W30 x D20 x H45 cm  
W11.8 x D7.9 x H17.7 in

PH-14 | 41430-000-

## Acrylic Phantom XAC



### Slab phantoms for radiation absorption and scattering measurement



SHOW MORE!



#### VARIATIONS


XAC-01	41430-000	30 x 30 x 0.1 cm/11.8 x 11.8 x 0.04 in
XAC-02	41431-000	30 x 30 x 0.2 cm/11.8 x 11.8 x 0.08 in
XAC-03	41432-000	30 x 30 x 0.3 cm/11.8 x 11.8 x 0.12 in
XAC-04	41433-000	30 x 30 x 0.4 cm/11.8 x 11.8 x 0.16 in
XAC-05	41434-000	30 x 30 x 0.5 cm/11.8 x 11.8 x 0.2 in
XAC-08	41435-000	30 x 30 x 0.8 cm/11.8 x 11.8 x 0.3 in
XAC-1	41436-000	30 x 30 x 1 cm/11.8 x 11.8 x 0.4 in
XAC-2	41437-000	30 x 30 x 2 cm/11.8 x 11.8 x 0.8 in
XAC-3	41438-000	30 x 30 x 3 cm/11.8 x 11.8 x 1.2 in
XAC-4	41439-000	30 x 30 x 4 cm/11.8 x 11.8 x 1.6 in
XAC-5	41440-000	30 x 30 x 5 cm/11.8 x 11.8 x 2.0 in
XAC-8	41441-000	30 x 30 x 8 cm/11.8 x 11.8 x 3.1 in
XAC-10	41442-000	30 x 30 x 10 cm/11.8 x 11.8 x 3.9 in

PH-16 | 41318-000,010 / 41319-000,010

## Contrast Detail Phantom



### Image evaluation in plain X-ray



SHOW MORE!



#### FEATURES

Four types of phantoms with different sizes and target types

- Hole 15: 41318-000      - Hole 10: 41318-010
- Rod 15: 41319-000      - Rod 10: 41319-010

#### DESCRIPTIONS

##### SIZE

- Hole 15: 15 x 15 holes of depth range from 1.0 to 8.0 mm (0.4 to 3.1 in)
- Rod 15: 15 x 15 rods of height range from 1.0 to 8.0 mm (0.4 to 3.1 in)
- Hole 10: 10 x 10 holes of depth range from 1.0 to 5.5 mm (0.4 to 2.2 in)
- Rod 10: 10 x 10 rods of height range from 1.0 to 5.5 mm (0.4 to 2.2 in)

##### SET INCLUDES

- 1 chart phantom
- storage case