

MW30

Head and Neck Palpation Trainer

New, innovative and ultra-realistic touch, designed to standardize palpation experiences


With realistic tactile material, this head and neck palpation trainer provides a sensation similar to human skin

Product Supervision

Junji Otaki,
professor and acting chief of the center,
Center for Medical Education,
Hokkaido University



KEY FEATURES



Swallowing

Swallowing motion to identify and assess the thyroid gland



Sites

Various deviations



Close-to-real-life sensation

Supple, resilient and delicate texture.

FEATURES


- 1 | Provides opportunities for hands-on training.
- 2 | Abnormal findings to complement training with healthy simulated patients
- 3 | Ensures that trainees share the same tactile sensation to specify the findings
- 4 | Maintain objective approaches for skills examination

SKILLS

- | Palpation of head and neck for lymph nodes, thyroid and salivary glands
- | Palpation skills can support clinical reasoning

CASES

	Sites	Findings
1	Posterior auricular nodes	Edema/Hard
2	Occipital nodes	N/A
3	Superficial cervical nodes	Edema/Palpable
4	Posterior cervical nodes	Edema/Palpable
5	Deep cervical nodes	Edema/Palpable
6	Pre-auricular nodes	Edema/Hard
7	Tonsillar nodes	Edema/Palpable
8	Submandibular nodes	Edema/Palpable
9	Submental nodes	Edema/Palpable
10	Supraclavicular nodes	Edema/Hard/Immovable
11	Subclavicular (infraclavicular) nodes	Edema/Hard/Immovable
12	Thyroid gland	Graves' disease
13	Submandibular gland	As hard as an enlarged lymph node
14	Parotid gland	As hard as an enlarged lymph node



MW31 Head and Neck Palpation Trainer for OSCE

The palpable nodes and glands in the standard type are distributed between two identical manikins for use in situations like student exams or skill checking.

Internal findings

- A** left side: 6/1/3/10 right side: 4/5/13/12
- B** left side: 3/13/ 9/5 right side: 4/11/5

DESCRIPTIONS

SPECIFICATIONS

Size: 29x29x47cm / 11.4x11.4x18.5 inch Weight: 6.5kg / 14.3lbs

SET INCLUDES

- 1 head and neck manikin with storage case
- 1 instruction manual

MATERIALS

Soft resin
Latex free