

41337-090

Breast Plate for Chest Phantom N-1

CT

X-ray

Simulates the absorption rate of breasts during female chest X-ray and demonstrates the influence of the nipple on image quality



SHOW MORE!



CHEST PHANTOMS

FEATURES

| This plate can be attached to PH-1 to study radiation absorption of the breasts and its effect on image quality.

APPLICATIONS

| CT
| Plain X-ray

DESCRIPTIONS

SET INCLUDES

1 breast plate

MATERIALS

urethane based resin (density: 1.06)

SPECIFICATIONS

Phantom size:
W47×D11×H55 cm
W18.5×D4.3×H21 in

Packing size:
W63×D50×H29 cm
W24.8×D19.7×H11.4 in

41337-080

Pneumonia Kit for N-1

CT

X-ray

Attachment for N-1 chest phantom to create true-to-life medical images of pneumonia



*The mediastinum and the pulmonary vessels are not included.

DESCRIPTIONS

SET INCLUDES

1 set of pneumonia attachments

MATERIALS

polyurethane (density:0.92)