

M66

# Medical Foot Care Simulator

## Simulate basic podiatric care

The maintenance of healthy feet makes a significant impact on the patient's quality of life. Training for simple trimming of nails and callosities as well as general education on basic foot care

### Product Supervision

Takako Egawa  
Dean, Professor, Kansai University of  
Nursing and Health Sciences  
Professor Emeritus, Kyoto University



Trimming and removal of the callosities



Trimming and removal of the corns

### FEATURES

- 1 | Realistic sensation and anatomically accurate
- 2 | Replaceable nails, corn and callosities

### SKILLS

- | Trimming and clipping of toe nails (thickened ingrown nail and thickened nail with ringworm)
- | Trimming and removal of callosities
- | Trimming and removal of corns
- | Foot assessment (demonstration)
- | Foot massage (demonstration)

## DESCRIPTIONS

### SPECIFICATIONS

Size: W21 x D21 x H28 cm / 8.3 x 8.3 x 11 inch Weight: 1.2g/2.9lbs

### SET INCLUDES

1	foot model with swivel stand	10	callosities
20	toe nails A (thickened ingrown nail)	10	corns
20	toe nails B (thickened nail with ringworm)	1	instruction manual

### REPLACEMENT PARTS

11273-020	1 replacement parts set for M66	11273-030 Replacement toe nail A	20 pcs
	-toe nail A	11273-040 Replacement toe nail B	20 pcs
	-toe nail B	11273-050 Replacement callosity	10 pcs
	-callosity	11273-060 Replacement corn	10 pcs
	-corn		10 pcs

\* Instruments are not included.

MW1



# Pitting Edema Simulator

## Five variations of pitting edema simulator

Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"

### Product Supervision

Junji Otaki  
Professor of Tokyo Medical University



Normal



1+ Light Pitting

2mm



2+ Moderate Pitting

4mm



3+ Deep Pitting

6mm



4+ Very Deep Pitting

8mm

### SKILLS

- | Differentiation of the pitting edema: four different levels of edemas and a normal case
- | Assessment of edemas on the lower legs

### FEATURES

- 1 | Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"
- 2 | The models can be attached to the lower leg of a simulated patients or a training manikin.
- 3 | 4 levels of pitting edema according to pressures and the depth of indentations.
- 4 | Repeatable and accurate, as if palpating actual skin surface to the fibula.
- 5 | Good combination with our nursing dolls

## DESCRIPTIONS

### SPECIFICATIONS

Size: W10 x D20 cm / W3.9 x D7.9 inch Weight: 0.15g/0.33lbs

### SET INCLUDES

5	variations of strap-on pitting edema model
1	instruction manual

### MATERIALS

Soft resin  
Latex free