

M66

Medical Foot Care Simulator

Simulate basic podiatric care

The maintenance of healthy feet makes a significant impact on the patient's quality of life. Training for simple trimming of nails and callosities as well as general education on basic foot care

Product Supervision

Takako Egawa
Dean, Professor, Kansai University of
Nursing and Health Sciences
Professor Emeritus, Kyoto University



Trimming and removal of the callosities



Trimming and removal of the corns

FEATURES

- 1 | Realistic sensation and anatomically accurate
- 2 | Replaceable nails, corn and callosities

SKILLS

- | Trimming and clipping of toe nails (thickened ingrown nail and thickened nail with ringworm)
- | Trimming and removal of callosities
- | Trimming and removal of corns
- | Foot assessment (demonstration)
- | Foot massage (demonstration)

DESCRIPTIONS

SPECIFICATIONS

Size: W21 x D21 x H28 cm / 8.3 x 8.3 x 11 inch Weight: 1.2g/2.9lbs

SET INCLUDES

1	foot model with swivel stand	10	callosities
20	toe nails A (thickened ingrown nail)	10	corns
20	toe nails B (thickened nail with ringworm)	1	instruction manual

REPLACEMENT PARTS

11273-020	1 replacement parts set for M66	11273-030 Replacement toe nail A 20 pcs
	-toe nail A 20 pcs	11273-040 Replacement toe nail B 20 pcs
	-toe nail B 20 pcs	11273-050 Replacement callosity 10 pcs
	-callosity 10 pcs	11273-060 Replacement corn 10 pcs
	-corn 10 pcs	

* Instruments are not included.

MW1



Pitting Edema Simulator

Five variations of pitting edema simulator

Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"

Product Supervision

Junji Otaki
Professor of Tokyo Medical University



Normal

1+ Light Pitting

2+ Moderate Pitting

3+ Deep Pitting

4+ Very Deep Pitting

2mm

4mm

6mm

8mm

SKILLS

- | Differentiation of the pitting edema: four different levels of edemas and a normal case
- | Assessment of edemas on the lower legs

FEATURES

- 1 | Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"
- 2 | The models can be attached to the lower leg of a simulated patients or a training manikin.
- 3 | 4 levels of pitting edema according to pressures and the depth of indentations.
- 4 | Repeatable and accurate, as if palpating actual skin surface to the fibula.
- 5 | Good combination with our nursing dolls

DESCRIPTIONS

SPECIFICATIONS

Size: W10 x D20 cm / W3.9 x D7.9 inch Weight: 0.15g/0.33lbs

SET INCLUDES

5	variations of strap-on pitting edema model
1	instruction manual

MATERIALS

Soft resin
Latex free