

Lumbar Puncture Simulator II

World's best training simulator for lumbar puncture and epidural anesthesia

Five types of puncture blocks enhance training with different types of challenges

Product Supervision

Takahiro Amano, M.D. Vice President Senior Vice-Dean, Postgraduate School Professor and Director, Center of Postgraduate Medical Education International University of Health and Welfare Honorary Director, Sanno Medical Center







FEATURES

- 1 | Anatomically accurate
- 2 | Realistic skin and tissue resistance during insertion
- 3 | Measure CSF fluid pressure, collect CSF fluid and practice epidural anesthesia procedures
- 4 | Simulates administration of epidural nerve block with realistic needle-tip sensation
- 5 | Simulates loss of resistance technique
- 6 | Model stand for either individual or group practice

SKILLS

- | Palpation of landmarks
- Lumbar puncture
- | CSF collection
- | CSF pressure measurement
- | Epidural anesthesia
- | Anatomical understanding



KEY FEATURES

- 5 types of puncture blocks for variable kinds
- | normal
- obese
- **I** senior
- I senior obese
- | epidural







DESCRIPTIONS

SPECIFICATIONS

Size: W33 x D21 x H30 cm W13 x D8.3 x H11.8 inch

SET INCLUDES 1 lumbar region model lumbar puncture blocks

> 2 normal CSF, 1 obesity CSF, 1 senior CSF 1 senior obesity CSF, 1 epidural

Weight: kg/ lbs

Power: -Power Consumption: -

1 lumbar region skin cover lumbar region support bases: sitting position, lateral position, base for team teaching 1 lumbar spine model

MATERIALS Soft resin Latex free

1 irrigation bag syringe

stand for irrigation bag instruction manual

1 storage case

RECOMMENDED DEVICES

REPLACEMENT PARTS

11348-090 1 LP block (normal CSF) 11348-110 1 LP block (obesity CSF)

11348-120 1 LP block (senior CSF)

11348-130 1 LP block (senior obesity CSF) 11348-140 1 epidural block 11348-150 1 lumbar region skin cover