

# Digital Rectal Examination Simulator



M92



Digital Rectal Examination Simulator is a unique trainer designed for training in palpation to identify rectal tumors and polyps.

## Features

1. 4 interchangeable rectum units for identification of lesion types.
2. Normal prostate or endocervix is also to be palpated.
3. Life-like rectum is palpable up to 7cm depth.
4. High-fidelity sphincter resistance and feeling of rectum wall allows excellent hands-on DRE training before examining actual patients.

## Set Includes

- 1 lumbar torso (unisex)
- 4 rectum units:
  - 1 normal
  - 2 rectal cancers
  - 1 rectum with polyps
- 1 prostate unit
- 1 endocervix unit
- 1 jar of vaseline
- 1 storage case
- 1 instruction manual

## Manikin Size

W37 x D28 x H30 cm


## Replacement Parts

- 11269-110
  - 11269-120
  - 11269-130
  - 11269-140
  - 11269-150
  - 11269-160
- 1 rectum unit #1 (normal)  
1 rectum unit #2 (cancer A)  
1 rectum unit #3 (cancer B)  
1 rectum unit #4 (rectum with polyps)  
1 prostate unit  
1 endocervix unit

# Prostate Examination Simulator



M53B



## Product Supervision

Dr. Osamu Yoshida, Professor Kyoto University Faculty of Medicine, Department of Urology  
Dr. Takashi Katayama, Professor Toyama Medical and Pharmaceutical University Department of Urology

Prostate Examination Simulator provides realistic training in prostate palpation skills with its soft life-like anus and rectum material.

## Features

1. Ten cases of prostate are set on a rotating table inside the body.
2. The cases are easily changeable.
3. Three patient positions are possible: supine, prone and lateral position.

## Cases

- 10 prostate cases:
- 2 normal
  - 2 enlargements
  - 5 carcinomas
  - 1 prostatitis

## Training Skills

- Digital palpation of the prostate

## Set Includes

- 1 male lumbar torso
- 1 rotating unit
- 10 prostate cases:
  - 2 normal, 2 enlargements,
  - 5 carcinomas and 1 prostatitis

## Manikin Size

W40 x D32 x H24 cm

## Replacement Parts

- 11272-040
- 1 replacement anus pad

- 1 lubricant gel
- 1 pair of gloves
- 1 storage case
- 1 instruction manual

rw-kyoto@kyotokagaku.co.jp