

MW18





Ultrasound-Guided PICC Training Simulator

New, innovative, human-like material for perfect visualization

A non-complicated design with replaceable puncture site and anatomically accurate features including the rib cage, sternum, bifurcation, plus visible tip placement in the superior vena cava.

Product Supervision

Joho Tokumine, M.D., Ph.D.,
Department of Anesthesiology,
Kyorin University School of Medicine



FEATURES

- 1 | Excellent image quality and visualization of the needle tip for ultrasound guided venous access
- 2 | Movable shoulder to demonstrate positioning
- 3 | Realistic flashback in needle provides confirmation for successful venous access
- 4 | Ribs and right clavicle provide anatomical understanding of correct PICC placement
- 5 | Anatomically correct bifurcation of the vein
- 6 | Simulation of cannula malposition


SKILLS

- | Ensure patient comfort and position
- | Hand hygiene and prepare sterile field
- | Identify puncture site
- | Ultrasound-guided insertion
- | Insertion of dilator and guide wire
- | Insertion of sheath
- | Insertion of catheter


KEY FEATURES



Ultrasound guided puncture



Movable shoulder



Malpositioning

DESCRIPTIONS

SPECIFICATIONS

Size: W40 x D15 x H60 cm Weight: 4.5kg/ 10lbs
W16 x D5.9 x H24 inch

SET INCLUDES

- | | |
|---|------------------------|
| 1 male upper torso with the right arm | 1 instruction manual |
| 2 PICC puncture pads | 1 storage case |
| 10 simulated blood (swab type) | |
| 1 jar | |
| 1 syringe | |

MATERIALS

Soft resin, hard resin
Latex free

REPLACEMENT PARTS

11398-010 2 PICC puncture pads

RECOMMENDED DEVICES

Ultrasound scanner and PICC instruments are not included.

