

MW32



NEW

NCPR Simulator

Manikin that facilitates Cardio-Pulmonary Resuscitation training of newborn baby



FEATURES

- 1 | Airway management training with observable bilateral/unilateral ventilation
- 2 | Soft and realistic umbilical cord capable of venous and arterial line placement, and blood sample collection
- 3 | Feeding tube placement confirmable through gastric fluid aspiration and
- 4 | Chest compression

SKILLS

Basic neonatal care

- | Changing body position
- | Hygiene
- | Palpation of the anterior fontanelle
- | Umbilical catheterization
- | Umbilical cord care
- | Tube feeding (nasal and oral) and gastric fluid aspiration
- | Clamping, tying and cutting of the umbilical cord
- | Suction (nasal, oral)

Emergency Simulation

- | Chest compression
- | Endotracheal intubation
- | BVM ventilation

*The product is not designed for bathing.

ANATOMY

- | ribs
- | bronchi
- | lungs
- | stomach
- | esophagus
- | umbilical cord

DESCRIPTIONS

SPECIFICATIONS

Size: 50 cm / 19.6 inch Weight: 2kg / 4.4 lbs

SET INCLUDES


- 1 newborn manikin
- 1 lubricant
- 1 instruction manual

MATERIALS



Soft resin
Latex free

REPLACEMENT PARTS

- 11400-010 1 mask for MW21, MW32
- 11400-020 1 chest skin for MW21, MW32
- 11400-030 4 umbilical cords for MW21, MW32
- 11229-050 1 lubricant



Palpation of the anterior fontanelle

Umbilical cord care
(one vein and two arteries)
umbilical catheterization and blood samplingTube Feeding (nasal and oral)
(tube placement confirmation with bubble sounds)

Airway suctioning (nasal and oral)



Chest compression

Endotracheal intubation
(confirmation of unilateral/bilateral ventilation)BVM ventilation
(rise of the chest)

RECOMMENDED DEVICES

Feeding tube: 5fr.
Suction catheter: 5fr.
Endotracheal tube: 3.0mm ID
Umbilical catheter: 5fr., 30cm
Laryngoscope: #0
BVM: neonatal

* Instruments are not included.

