

M43B/M43B-B



Lumbar Puncture Simulator IIA


World's best training simulator for lumbar puncture and epidural anesthesia

Five types of puncture blocks enhance training with different types of challenges



Product Supervision

Takahiro Amano, M.D.
Vice President
Senior Vice-Dean, Postgraduate School
Professor and Director, Center of Postgraduate
Medical Education
International University of Health and Welfare
Honorary Director, Sanno Medical Center



FEATURES

- 1 | Anatomically accurate
- 2 | Realistic skin and tissue resistance during insertion
- 3 | Measure CSF fluid pressure, collect CSF fluid and practice epidural anesthesia procedures
- 4 | Simulates administration of epidural nerve block with realistic needle-tip sensation
- 5 | Simulates loss of resistance technique
- 6 | Model stand for either individual or group practice

SKILLS

- | Palpation of landmarks
- | Lumbar puncture
- | CSF collection
- | CSF pressure measurement
- | Epidural anesthesia
- | Anatomical understanding

KEY FEATURES

5 types of puncture blocks for variable kinds

- | normal
- | obese
- | senior
- | senior obese
- | epidural



DESCRIPTIONS

For product codes of the alternate skin tone versions of this model, please refer to www.kyotokagaku.com.

SPECIFICATIONS

Size: W33 x D22 x H30 cm Weight: 3.15kg / 6.9lbs
W13 x D8.7 x H11.8 inch

MATERIALS

Soft resin
Latex free

RECOMMENDED DEVICES

-

SET INCLUDES

	lumbar puncture blocks	M43B	M43B-B
1 lumbar region model	2 normal CSF	✓	✓
1 lumbar region skin cover	1 obesity CSF	✓	✓
1 two way support base	1 senior CSF	✓	✓
1 base for team teaching	1 senior obesity CSF	✓	✓
	1 epidural	✓	✓
	1 lumbar spine model	✓	✓

- 1 irrigation bag
- 1 syringe
- 1 stand for irrigation bag
- 1 instruction manual
- 1 storage case

REPLACEMENT PARTS

11348-090	1 LP block (normal CSF)
11348-110	1 LP block (obesity CSF)
11348-120	1 LP block (senior CSF)
11348-130	1 LP block (senior obesity CSF)
11348-140	1 epidural block
11348-150	1 lumbar region skin cover
11348-160	1 supporting stand