

M43E

Ultrasound Compatible Lumbar Puncture/ Epidural Simulator

Lumbar puncture training with ultrasound guidance

Anatomically correct puncture block offers excellent ultrasound images






FEATURES

- 1 | Durable and replaceable parts for repetitive practice
- 2 | Adjustable to an upright or lateral position.
- 3 | Translucent block makes the needle trace visible.
- 4 | Interchangeable parts for other landmark and fluoroscopic procedure training

SKILLS

- | Ultrasound-guided lumbar puncture
- | Ultrasound-guided epidural anesthesia
- | CSF collection and CSF pressure measurement

IMAGES



DESCRIPTIONS

SPECIFICATIONS

Size: W33 x D21 x H30 cm
W13 x D8.3 x H11.8 inch

SET INCLUDES

- | | |
|---|----------------------|
| 1 lumbar region model | 1 irrigation bag |
| 1 ultrasound lumbar puncture/epidural block | 1 syringe |
| 1 ultrasound lumbar region skin cover | 1 instruction manual |
| 1 two way support base | |

MATERIALS

Soft resin, polyurethane elastomer
Latex free

REPLACEMENT PARTS

- | | |
|-----------|--|
| 11348-190 | 1 ultrasound lumbar puncture /epidural block |
| 11348-230 | 1 ultrasound lumbar region skin cover for M43E |
| 11348-160 | 1 supporting stand |

RECOMMENDED DEVICES

Ultrasound machine and needles are not included.