

M53B



Prostate Examination Simulator

Prostate Examination Simulator provides realistic training in prostate palpation skills with its soft life-like anus and rectum material.



FEATURES

- 1 | Ten cases of prostate are set on a rotating table inside the body
- 2 | The cases are easily changeable
- 3 | Three patient positions are possible: supine, prone and lateral position

SKILLS

Digital palpation of the prostate

Cases

2 normal, 2 enlargements, 5 carcinomas, 1 prostatitis



DESCRIPTIONS

SPECIFICATIONS

Size: W30 x D36 x H24 cm
W11.8 x D14.2 x H9.4 inch

Weight: 5.1kg/ 10.9lbs

MATERIALS

Soft resin, Hard resin
Latex free

SET INCLUDES

- 1 male lumbar torso
- 1 rotating unit
- 10 prostate cases: 2 normal, 2 enlargements, 5 carcinomas and 1 prostatitis
- 1 lubricant gel
- 1 pair of gloves
- 1 storage case
- 1 instruction manual