

M56B



# Pediatric Patient Care Simulator

## A comprehensive day-to-day, long-term care manikin for pediatric patients

Provides a variety of features for practicing essential skills and simulates pediatric patient care scenarios for nurses and caregivers

### Product Supervision

Medical Association of Osaka Prefecture  
Pediatric Medical Care Exploratory Committee



### FEATURES

- 1 | Life-like 6 years old model
- 2 | Able to sit in an upright position without support.
- 3 | Oral, nasal and tracheal airway suction interchangeable female and male genitalia for catheterization training.

### SKILLS

- | Positioning
- | Patient handling
- | Airway suction:
  - oral, nasal or through a tracheostomy opening
- | Percutaneous Endoscopic Gastrostomy (PEG) care
- | Feeding tube management
- | Male & female catheterization

### IMAGES



Airway suction



Airway suction




Airway suction



Feeding tube management



PEG care



Positioning, patient handling



Male &amp; female catheterization



Anatomical understanding

### DESCRIPTIONS

#### SPECIFICATIONS

Size: 125cm / 49.2 inch Weight: 13kg / 28.6lbs

#### SET INCLUDES

- |                             |                                    |
|-----------------------------|------------------------------------|
| 1 pediatric manikin         | 1 strap-on IV venipuncture trainer |
| 1 head & neck anatomy plate | 1 uniform (gym clothes)            |
| 1 male & female genitalia   | 1 set of catheters (sample)        |
| 1 syringe                   |                                    |
| 1 irrigation bag            |                                    |

#### MATERIALS

Soft resin, hard resin  
Latex free

#### RECOMMENDED DEVICES

-

#### REPLACEMENT PARTS

- 11355-020 head and mask explanation panel
- 11355-030 1 valve for catheterization
- 11355-040 replacement mouth/ nose part

\* Instruments are not included.

