M57B



Suture Evaluation Simulator

Score your suturing

Evaluate simple interrupted skin suturing techniques by quantitative assessment criteria. Objective evaluation score and record of learning curve would encourage and motivate trainees to improve skills



Introductory Ultrasound Training Block "REAL VESSEL"

Provides training in hands-eye coorditaion and basic skills in ultrasound-guided venous access.



FEATURES

- 1 | 2 simulated vessel lines: straight and curving.
- 2 | Lines have slope represent vessels with different depth.
- 3 | Vessel wall yields under pressure of a needle tip.

SKILLS

- | Visualization and localization of the vessels.
- Probe manipulation.
- | Basics for ultrasound-guided vessel access.

SPECIFICATIONS Size: W10.5 x D26 x H5.5 cm W4.2 x D10.3 x H2.2 in.

SET INCLUDES

- 5 vaccination training pads
- 5 pad holders
- 5 strap-on belts
 - instruction manual

DESCRIPTIONS -

MATERIALS Soft resin Latex free

RECOMMENDED DEVICES

REPLACEMENT PARTS

11454-010MW64 Pad for Vaccination Training Set (Set of 5)11454-020MW64 Pad for Vaccination Training Set

KYOTO KAGAKU