



M63 Female

M59 Male

Multipurpose newborn for skills training.

This training model is a life-size seamless manikin of a 0-4-week-old infant, which has the weight and balance of the newborn as well as palpable anterior and posterior fontanel.

Features

1. Life-like, life-size newborn model
2. Non-rigid neck flexibility
3. The fonticulus anterior and posterior are palpable
4. The limbs can bend and stretch.
5. Waterproof: the model has a seamless, one-piece structure
6. The umbilical cord is removable.

Training Skills

- Bathing
- Umbilical cord care
- Bottle-feeding
- Ear and nose care
- Diaper and clothing changing
- Physical measurement



Specifications :

Set Includes:	Manikin Size:
1 newborn manikin	50cm
1 umbilical cord	
1 feeding bottle	
1 baby clothes	
1 cloth cover	

* A baby bottle is not included.



Life-size newborn manikin for basic vital sign assessments

Newborn Vital Signs Examination Simulator is an effective tool for nurses to learn how to take vital signs of a newborn child.

Features

1. Arterial pulse in wrists, neck and temples
2. Rectal temperature can be measured.
3. Heart sounds are recorded from actual patients and synchronized with arterial pulses.
4. All vital signs can be controlled with the control unit.

Training Skills

- Vital signs examination
 - Heart sounds
 - Respiration
 - Pulses:
 - superficial temporal, carotid and radial arteries (bilateral)
 - Body temperature (rectal thermometry)

Product Supervision

Keiko Takahashi, professor
College of Nursing Science, Fukui Prefecture University



Palpation of Pulses



Auscultation

Specifications :

Set Includes:	Manikin Size:
1 newborn manikin	50cm
1 control box	
1 thermometer	
1 baby clothes	
1 storage case	

* A stethoscope is not included.