

M72

Bathing and In-patient Care Simulator (Male)

Elderly whole-body manikin for basic patient care

This manikin facilitates training for various procedures essential for elderly patient care, such as bathing, hair care, positioning and oral care.



FEATURES

- 1 | Flexible and seamless joints
- 2 | Waterproof body
- 3 | Allows practice of shampooing, brushing and drying patient's hair
- 4 | Skill training for such patient handling as positioning, carrying or transferring on and off of bed



SKILLS

Positioning
Patient handling
Range of joint-motion training
Hair care
Oral care
Diaper changing, clothing
Bathing



DESCRIPTIONS

SPECIFICATIONS

Size: 168 cm / 66 in Weight: 14 kg / 30 lbs

SET INCLUDES

- 1 manikin
- 1 denture
- 1 wig
- 1 underpants

MATERIALS

Soft resin
Latex free

REPLACEMENT PARTS

11274-020 wig for M72