

# **M73**

# **Baby Patient Care Model**

7-10 months old

---

## **Instruction Manual**

## Index

This training model is a life-size seamless manikin of a 7-10-month old child.

The simulator is designed for training in nursing care.

Do not use for improper purpose.

Please read the instruction carefully before using

## Specifications

1. Life-like, life-size 7-10 months old model
2. The model is designed to provide the feeling of approx.9kg weight for patient handling practice.  
(9 kg is not the actual weight)

### Set includes

1 manikin

1 set of pajamas

## DOs and DON'Ts

### DOs

#### Handle with care

The model is made of special composition of soft resin. Please handle it with the utmost care like an actual baby at all times. Store the manikin with its clothes on, covered by the supplied cloth, away from heat, moisture and direct sunlight.

### DON'Ts

The model is not designed for bathing. **NEVER immerse the manikin in water.**

Don't mark on the model with pen or leave any printed materials in contacted with their surface.

**Ink marks on the models will not be removable.**

## KYOTO KAGAKU CONTACT & ORDERING INFORMATION

<http://www.kyotokagaku.com> e-mail: [rw-kyoto@kyotokagaku.co.jp](mailto:rw-kyoto@kyotokagaku.co.jp)

#### ■ Main Office and Factory

15 Kitanechoya-cho Fushimi-ku Kyoto 612-8388, Japan  
Telephone : 81-75-605-2510  
Facsimile : 81-75-605-2519

#### ■ KyotoKagaku America Inc.

USA,Canada,and South America  
3109 Lomita Boulevard, Torrance, CA 90505 ,USA  
Telephone : 1-310-325-8860  
Facsimile : 1-310-325-8867