

M74



# Gluteus Maximus Injection Simulator

## IM training for the three ways of injection

This simulator is designed to practice for the three ways of injection into the hip and the buttocks. The incorporated indicator shows the accuracy of needle placement and its depth

### Product Supervision

Dr. Hiroshi Koyama,  
Medical Director, General Internal Medicine  
National Hospital Organization Kyoto Medical  
Center



### FEATURES

- 1 | The life-sized lower torso manikin consist of three layers: skin, muscular and skeleton
- 2 | Greater trochanter, anterior superior iliac spine and scrum can be palpated to identify the injection sites
- 3 | Three methods for intramuscular injection: Quarter-Trichotomy, Hochstetter Method and Clark Method
- 4 | The indicator shows the accuracy of the needle placement and its depth
- 5 | Correct needle placement can be confirmed by being able to inject liquid
- 6 | No needle trace remains

### SKILLS

- | Intramuscular injection



Hochstetter Method



Clark / Sachtleben Method

### DESCRIPTIONS

#### SPECIFICATIONS

Size: W32 x D35 x H20 cm      Weight: 3.2kg / 7.05lbs  
W12.6 x D14 x H8 in.

#### SET INCLUDES

- 1 | lower torso manikin
- 1 | stand for Hochstetter method
- 1 | indicator
- 1 | set of sensor pads
- 1 | instruction manual

#### MATERIALS

Soft resin  
Latex free

#### RECOMMENDED DEVICES

IM needle 21-23G

#### REPLACEMENT PARTS

- 11213-010 skin
- 11213-020 sensor for puncture pad
- 11213-030 sensor for Hochstetter's site
- 11213-040 sensor for Clark's site