


NEW

Airway Suction Trainer II

Life-like anatomy and realistic tactile feedback for intermittent airway suction training.



Product No. M85A
Part No. 11229-300



FEATURES

- | 01 Transparent anatomical structure to monitor the catheter
- | 02 Hands-on suction training using simulated sputum.
- | 03 The internal anatomy of the nasal cavity is faithfully reproduced, allowing learners to understand the key points of catheter insertion including the insertion length

SKILLS

Intermittent airway suction (oral, nasal, endotracheal)
Suction training on patients with a bite block or an inserted endotracheal tube

[Recommended devices]

Suction catheter 14Fr

Tracheostomy tube

Oral cavity / Tongue / Oral and nasal pharynx / Epiglottis /
Trachea / Esophagus

Redesigned for easier disassembly and cleaning



DESCRIPTIONS

SPECIFICATIONS

Adult head and chest manikin
Dimensions: W41×D28×H21cm
Weight: 1.4kg

Materials: Soft resin

SET INCLUDES

1	head and neck manikin with airway anatomy	1	simulated sputum
2	sample catheter 14Fr	1	lubricant
1	syringe 50mL	1	storage case

RELATED PRODUCTS

MW8A Airway Suction Simulator Plus
MW5 11229-200 Tube Feeding Panel
MY-102 12746-020 Tracheotomy tube II
11229-060 Oral cover for suctioning



KYOTO KAGAKU



rw-kyoto@kyotokagaku.co.jp
www.kyotokagaku.com

July, 2026