

M92



Digital Rectal Examination Simulator

Digital Rectal Examination Simulator is a unique trainer designed for training in palpation to identify rectal tumors and polyps.



FEATURES

- 1 | 4 interchangeable rectum units for identification of lesion types.
- 2 | Normal prostate or endocervix is also to be palpated.
- 3 | Life-like rectum is palpable up to 7cm depth
- 4 | High-fidelity sphincter resistance and feeling of rectum wall allows excellent hands-on DRE training before examining actual patients.

SKILLS

Digital palpation of the prostate

DESCRIPTIONS

SPECIFICATIONS

Size: W37 x D28 x H30 cm

REPLACEMENT PARTS

11269-210 1 rectum unit #1 (normal)
 11269-220 1 rectum unit #2 (cancer A)
 11269-230 1 rectum unit #3 (cancer B)
 11269-240 1 rectum unit #4
 (rectum with polyps)
 11269-150 1 prostate unit
 11269-160 1 endocervix unit

SET INCLUDES

1 lumbar torso (unisex)
 4 rectum units:
 1 normal
 2 rectal cancers
 1 rectum with polyps

1 prostate unit
 1 endocervix unit
 1 jar of Vaseline
 1 storage case
 1 instruction manual

M53B



Prostate Examination Simulator

Prostate Examination Simulator provides realistic training in prostate palpation skills with its soft life-like anus and rectum material.



FEATURES

- 1 | Ten cases of prostate are set on a rotating table inside the body
- 2 | The cases are easily changeable
- 3 | Three patient positions are possible: supine, prone and lateral position

SKILLS

Digital palpation of the prostate

Cases

2 normal, 2 enlargements, 5 carcinomas, 1 prostatitis

DESCRIPTIONS

SET INCLUDES

1 male lumbar torso
 1 rotating unit
 10 prostate cases: 2 normal,
 2 enlargements, 5 carcinomas
 and 1 prostatitis
 1 lubricant gel
 1 pair of gloves
 1 storage case
 1 instruction manual

REPLACEMENT PARTS

11272-040 1 anus pad

SPECIFICATIONS

Size: W30 x D36 x H24 cm
 W11.8 x D14.2 x H9.4 inch

Weight: 5.1kg/ 10.9lbs

MATERIALS

Soft resin, Hard resin
 Latex free