

M93UB/U

11347-300,500



CVC Insertion Simulator II

Conduct a variety of CVC insertion and management simulations with the assistance of 3 interchangeable pads




Product Supervision

Masahiro Tanabe, M.D., Ph.D., Professor Director General Medical Education Center Chiba University School of Medicine

Kinya Sando, M.D., Ph.D., Professor Director Department of Human Dietetics Graduate School of Human Science Osaka Shoin Women's University

Masanori Hoki, M.D., Ph.D., Professor Department of Human Dietetics Graduate School of Human Science Osaka Shoin Women's University

Joho Tokumine, M.D.Ph.D., Department of Anesthesia, Kawatetsu Chiba Hospital



FEATURES

- 1 Realistic ultrasound-guided CVC
- 2 Landmark puncture pad with anatomically correct vein bifurcations simulates mechanical complications including pneumothorax, mislodging and artery puncture*
- 3 Transparent anatomical block for anatomical understanding and guide wire manipulation
- 4 Introductory ultrasound training block to acquire basics of ultrasound-guided venous access
- 5 Both internal jugular and subclavian (axillary) veins are accessible

SKILLS

- | Ultrasound-guided CVC
- | Landmark-guided CVC*
- | Ultrasound-guided venous access
- | Prevention of mechanical complications

*Landmark pad is included only in M93UB set

KEY FEATURES



Transparent anatomical block



Ultrasound training block



Ultrasound puncture pad

DESCRIPTIONS

SPECIFICATIONS

Size: W40 x D20 x H 34 cm
W15.8 x D7.9 x H 13.4 in.
Weight: 2.3kg/5lbs

MATERIALS

Soft resin
Latex free

SET VARIATIONS

COMMON

- 1 male upper torso manikin
- 1 skin for cannulation training
- 1 transparent anatomical block
- 1 introductory ultrasound training block
- 1 irrigation bottle
- 1 syringes (2.5ml, 50ml)
- 1 sample needle
- 1 instruction manual

M93UB

- 1 landmark puncture pad
- 1 ultrasound puncture pad
- 10 cotton swab red
- 10 cotton swab blue
- 2 plastic jars
- 1 air bulb

M93U

- 2 ultrasound puncture pad
- 1 plastic jar

REPLACEMENT PARTS

- 11347-280 2 landmark puncture pads
- 11347-170 2 ultrasound puncture pads
- 11347-210 2 introductory ultrasound training blocks
- 11347-320 irrigation bottle
- 11388-400 10 cotton swab red
- 11402-010 10 cotton swab blue

Ultrasound scanner and needles are not included.