MW50

Strap-on Venipuncture Simulator

Cost-effective training unit for venipuncture

Strap-on Venipuncture Simulator is a simple, economical training kit for intravenous injection. This trainer is to be strapped onto the upper arm of an SP (Standardized Patient) or a life-size training manikin.



FEATURES

- 1 | Realistic tissue and vein wall resistance
- 2 | Vein tube can be filled and refilled easily using a soft plastic bottle that works as a manual pump.

RECOMMENDED DEVICES

 3 Three vein tubes are embedded: 2 normal and 1 thin.

SKILLS

| Blood collection and intravenous injection

DESCRIPTIONS

ISPECIFICATIONS Size: W25 x D11 x H3 cm W9.8 x D4.3 x H1.2 inch

SET INCLUDES

- 5 puncture pads (vein tube embedded)
- 5 supporting frames
- 5 belts
- 5 injection bottles
- 5 syringes

instruction manual

MATERIALS Soft resin Latex free

REPLACEMENT PARTS

11431-010 5 replacement puncture pads

OPTIONAL PARTS

11430-010 Arm Vessels Anatomy Sleevelet 11388-400 Cotton Swab for Coloring (Red) 10 pcs

M94

Intradermal Injection Training Simulator

Wheal can be shown by correct injection

This simulator provides training for intradermal injection to the inner forearm. Swelling of the skin will form when proper injection procedures are followed.

FEATURES

- 1 | Life-like touch and resistance of the injection site resembles the
- 2 | Injection pads are easily changeable.
- 3 | No needle trace remains.





SKILLS

| Intradermal Injection

ISPECIFICATIONS

Size: W9 x D20 x H8 cm W3.5 x D7.9 x H3.1 inch

Weight: 0.3 kg/ 0.7 lbs

SET INCLUDES

- 5 puncture pads
- forearm model
- 1 supporter stand storage case
- 1 instruction manual

*Needle not included.

Needle 29G

RECOMMENDED DEVICES

DESCRIPTIONS MATERIALS

> Soft resin Latex free

REPLACEMENT PARTS

11346-010 1 replacement puncture pad