

M95

11350-000

Simulator Pediatric IV Hand



Set of pediatric hands for practicing IV

Become confident in safe insertion of intravenous cannulas on small infant hands with thin and soft veins

FEATURES

- 1 | Fingers and dorsum of hand are flexible, able to be gripped, and repositionable as in real-life procedures
- 2 | Two puncture sites on the back of the hand
- 3 | Practice dressing and securing the cannula

SKILLS

- | Pediatric dorsal veins injection
- | Blood collection
- | Intravenous infusion

SPECIFICATIONS

Size: W18 x D30 x H9.5 cm / W7.1 x D11.8 x H3.7 in. (one-year-old)
W18 x D36 x H10 cm / W7.1 x D14.2 x H3.9 in. (three-year-old)

Weight: 0.6 kg / 1.3 lbs (one-year-old)
0.7 kg / 1.5 lbs (three-year-old)

Power: AC100-240V
Power Consumption: 5W

MATERIALS

Soft resin, hard resin
Latex free

SET INCLUDES

- | | | | |
|----|--------------------------------------|---|--------------------|
| 2 | injection arm-hands (one-year-old) | 1 | spoon |
| 2 | injection arm-hands (three-year-old) | 1 | plastic jar |
| 2 | stands | 1 | plastic beaker |
| 1 | circulation pump | 1 | storage case |
| 10 | cotton swab red | 1 | instruction manual |

DESCRIPTIONS

REPLACEMENT PARTS

- | | |
|-----------|---|
| 11350-010 | 1 replacement injection arm-hand (three-year-old) |
| 11350-020 | 1 replacement injection arm-hand (one-year-old) |
| 11388-400 | cotton swab red |

RECOMMENDED DEVICES

indwelling needle: 23G
butterfly needle: 23G

* Needles are not included.

M43D



11348-400,-M,-D


Pediatric Lumbar Puncture Simulator II

Acquire skills for safe pediatric lumbar puncture

The full body is soft and bendable to facilitate realistic body positioning required for a safe and successful lumbar puncture



Level of posterior superior iliac crests



FEATURES

- 1 | Soft, seamless baby manikin allows practice in safe patient positioning and holding
- 2 | Two types of puncture pads: deep or shallow
- 3 | Realistic needle tip resistance
- 4 | Anatomically correct landmarks for palpation

SKILLS

- | Patient positioning
- | Palpation of landmarks
- | Lumbar puncture and CSF collection
- | CSF pressure measurement

Product Supervision

Takahiro Amano, M.D.
Vice President
Senior Vice-Dean, Postgraduate School
Professor and Director, Center of Postgraduate
Medical Education International University of Health
and Welfare,
Honorary Director, Sanno Medical Center
Gregory A. Plotnikoff, M.D., M.T.S., F.A.C.P.
Medical Director, Abbott Northwestern Hospital

SPECIFICATIONS

Size: W43 x D34 x H20 cm
W17 x D13.4 x H7.8 in.

Weight: 1.5 kg / 3.3 lbs

- | | | | |
|---|------------------------------|---|--------------------|
| 1 | baby manikin (7-10 months) | 1 | syringe (50mL) |
| 2 | pediatric LP block (shallow) | 1 | instruction manual |
| 2 | pediatric LP block (deep) | | |
| 1 | irrigation bag | | |
| 1 | stand for irrigation bag | | |

DESCRIPTIONS

MATERIALS

Soft resin
Latex free

RECOMMENDED DEVICES

22G-24G spinal needle

* Needles are not included.

REPLACEMENT PARTS

- | | | | |
|--------------|---------------------------------------|-------------|------------------------------------|
| 11348-100-04 | stand for irrigation bag | 11348-420 | 2 pediatric LP block light (deep) |
| 11348-100-05 | irrigation bag | 11348-420-M | 2 pediatric LP block medium (deep) |
| 11348-410 | 2 pediatric LP block light (shallow) | | |
| 11348-410-M | 2 pediatric LP block medium (shallow) | 11348-420-D | 2 pediatric LP block dark (deep) |
| 11348-410-D | 2 pediatric LP block dark (shallow) | | |