

Ultrasound Guided Arterial Puncture Training Pads

Enhanced Training with the M99: Master Ultrasound-Guided Radial Arterial Puncture

Boost your confidence and skill in arterial blood sampling



Ultrasound Puncture pads



Training with M99 Arterial Puncture Training Wrist


FEATURES

- 1 | Palpable pulsation of arterial pulse
- 2 | Realistic needle resistance against tissue and artery wall
- 3 | Natural counter-flow
- 4 | Radial arterial line placement can be simulated
- 5 | No puncture trace remains on the injection site
- 6 | Leak-free connections


SKILLS

- | Ultrasound-guided
 - radial artery puncture
 - artery blood sampling
 - radial line placement


IMAGES



Transverse Plane



Longitudinal Plane



IABP Measurement
(Arterial pressure values are not simulated)

DESCRIPTIONS

REPLACEMENT PARTS

- 11351-110 2 ultrasound puncture pads
- 11351-120 2 ultrasound puncture inner pads
- 11351-130 2 ultrasound puncture outer pads
- 11388-400 10 cotton swabs red

MATERIALS

- Soft resin
- Latex free

RECOMMENDED DEVICES

- AGB Kit 22-25G

Next page : M99 Arterial Puncture Training Wrist

Ultrasound scanner and needles are not included.