

MW1



Pitting Edema Simulator

Five variations of pitting edema simulator including a normal one

Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"

Product Supervision

Junji Otaki
Professor of Tokyo Medical University



FEATURES

- 1 | Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"
- 2 | The models can be attached to the lower leg of a simulated patients or a training manikin.
- 3 | 4 levels of pitting edema according to pressures and the depth of indentations.
- 4 | Repeatable and accurate, as if palpating actual skin surface to the fibula.

SKILLS

- | Differentiation of the pitting edema: four different levels of edemas and a normal case
- | Assessment of edemas on the lower legs

KEY FEATURES



Normal

1+ Light Pitting

2+ Moderate Pitting

3+ Deep Pitting

4+ Very Deep Pitting

DESCRIPTIONS

SPECIFICATIONS

Size: 10x23cm / 3.9x9 inch

SET INCLUDES

5 variations of strap-on pitting edema model 1 piece each

MATERIALS

Soft resin
Latex free