


MW20

Prostate and Rectal Simulator

2-in-1 simulator for prostate and rectal exams

The Prostate and Rectal Simulator includes 9 types of prostates and 4 types of rectums which can easily be changed. The lower torso manikin can be positioned in 3 different ways.



FEATURES

- 1 | Nine types of prostates simulate different scenarios
- 2 | Four rectum units
- 3 | Training with multiple prostates and rectal cases are concurrent and time efficient
- 4 | Lateral, supine and prone positioning
- 5 | 3 revolving units with 9 prostates can easily be changed.

SKILLS

- | Digital-examination of the prostate and rectum
- | Insertion and the use of the anal-scope and the proctoscope

KEY FEATURES

Cases: Nine types of Prostates



- | | | |
|------------------|----------------|----------------|
| #1 normal | #4 prostatitis | #7 carcinoma 3 |
| #2 enlargement 1 | #5 carcinoma 1 | #8 carcinoma 4 |
| #3 enlargement 2 | #6 carcinoma 2 | #9 carcinoma 5 |

Cases: Four Rectum Units



- | | |
|--------------------|--|
| #1 normal | #3 large carcinoma |
| #2 small carcinoma | #4 polyp combines with small carcinoma |

DESCRIPTIONS

SPECIFICATIONS

Size: 34 x 33 x 24 cm/13.4 x 13 x 9.4 inch

SET INCLUDES

1 lower torso manikin	4 rectum units
3 prostate units	1 Vaseline
(include 9 types of prostates)	1 storage case

MATERIALS

Soft resin
Latex free

REPLACEMENT PARTS

11399-010 prostate unit A (normal, enlargement 1&2)
11399-020 prostate unit B (prostatitis, carcinoma 1&2)
11399-030 prostate unit C (carcinoma 3,4&5)
11399-040 rectum (normal)
11399-050 rectum (small carcinoma)
11399-060 rectum (large carcinoma)
11399-070 rectum (polyp and small carcinoma)
11399-080 skin