


MW21B



NEW

NCPR Simulator Plus



Intravenous catheter insertion (vein localization via transillumination is possible)

Added functions of PICC and IVC insertion to the standard NCPR Simulator

FEATURES

- 1 | Transillumination of the right hand and foot for IVC insertion
- 2 | Accessible veins on the right arm and leg for PICC insertion

SKILLS

In addition to NCPR Simulator, NCPR Simulator Plus trains: IV access with transilluminator guidance

Accessible veins:

- | right arm basilic vein
- | right hand dorsal vein
- | right leg saphenous vein
- | right leg popliteal vein

*The product is not designed for bathing.

SPECIFICATIONS

Size: 50 cm / 19.6 inch
Weight: 2.6kg / 5.7 lbs

SET INCLUDES

- 1 newborn manikin
- 1 lubricant
- 1 instruction manual

MATERIALS

Soft resin
Latex free

DESCRIPTIONS

REPLACEMENT PARTS

- 11400-010 mask
- 11400-020 chest skin
- 11400-030 umbilical cord
- 11400-040 1 right arm skin
- 11400-050 1 right leg skin
- 11400-310 4 vein tubes for right arm
- 11400-320 4 vein tubes for right leg
- 11400-080 abdomen
- 11388-400 10 simulated blood sticks (swab type)
- 11229-050 lubricant

RECOMMENDED DEVICES

- Feeding tube: 5 fr.
- Suction catheter: 5 fr.
- Endotracheal tube: 3.0mm ID
- Peripheral venous catheter: 28G x 30 cm
- Umbilical catheter: 5fr., 30cm
- Laryngoscope: #0
- BVM: neonatal

* Instruments are not included.

MW33



NEW

ELBW Infant Simulator



26-29 weeks gestation



Extremely low birth weight infant model offers the opportunities to practice fundamental skills in NICU

FEATURES

- 1 | A male infant of 26-29 weeks gestational, 900g weight
- 2 | Resuscitation training including chest compression and endotracheal intubation
- 3 | Life-like softness and delicateness
- 4 | Articulated head and neck

SKILLS

- | Changing body position, hygiene care in an incubator
- | Whole body assessment and measurement
- | Oral gastric and Nasogastric Tube placement
- | Insertion of urinary catheter
- | Bag valve mask resuscitation
- | Chest compression
- | Endotracheal intubation
- | Nasal CPAP

ANATOMY

- | bronchi
- | esophagus
- | stomach
- | articulated head and neck

*The product is not designed for bathing.

SPECIFICATIONS

Size: 35cm / 13.7 inch
Head circumference: 28cm / 11 inch

SET INCLUDES

- 1 ELBW manikin
- 1 lubricant
- 1 instruction manual

Weight: 0.9kg / 2 lbs

MATERIALS

Soft resin
Latex free

REPLACEMENT PARTS

- 11413-010 1 mask for MW33
- 11229-050 1 lubricant

RECOMMENDED DEVICES

- Feeding tube: 3fr.
- Suction catheter: 3fr.
- Endotracheal tube: 2.0-2.5 mm ID
- Laryngoscope: #00
- Urine catheter: 6fr.
- BVM: neonatal (small)

* Instruments are not included.

