


MW23

Adult Diaper Changing Trainer

Promotes clear visual understanding in adjusting a diaper in order to ensure patient's comfort and prevent leakage.

Appropriate use of a diaper improves patient's Activities of Daily Living (ADL) while preventing bedsores and their aggravation




FEATURES


- 1 | Landmarks for a diaper adjustment.
- 2 | Simulated urine or watery feces can be used.
- 3 | Transparent lower torso allows for direct observation of internal view of the diaper.
- 4 | Both male and female patients can be simulated.
- 5 | Lightweight and easy-to handle.

SKILLS

- | Elimination assistance (Diaper-Changing)



male genitalia



DESCRIPTIONS

SPECIFICATIONS

Life-sized lower torso: W31 x D50 x H21cm, 0.5kg

SET INCLUDES

1 lower torso manikin	1 simulated watery feces
1 male genitalia	1 bottle
1 irrigation bag	1 manual pump
1 connection tube	1 set of simulated urine
1 syringe	(10 swab set)
1 beaker	

REPLACEMENT PARTS

11381-800 simulated urine (10 swab set)
 11229-070 simulated watery feces (set of 3 packs)
 11402-010 10 cotton swab blue