

MW26





Basic Patient Care Simulator "Keiko"

Meet Keiko :

Patient simulator of geriatric care for trainees including home-care nurses and caregivers

Keiko simulates patients under long-term cares in a variety of health care settings



Feeding

- | Tube feeding (oral/nasal)

Exercise

- | Body positioning
- | Patient transfer

Wound Care

- | Pressure sore care

Pain Management

- | Hot and Cold Compressions

Elimination

- | Enema
- | Urinary catheterization
- | Stoma care
- | Bedpan practice

Hygiene

- | Oral care
- | Bed bathing
- | Change of clothes
- | Perineal care
- | Hair care*
- | Grooming and appearance

Respiratory Care

- | Oxygen inhalation
- | Airway suction
- | Postural drainage
- | Tracheotomy care

Others

- | Postmortem care

Drug Administration

- | Suppository

* Optional



KEY FEATURES

Basic Patient Care

Oral Care

- Oral care on the bed and wheel chair
- Special care for patients with feeding tubes



- | Oral assessment
- | Teeth brushing
- | Denture care

Pressure sores care

- Learn stages of pressure sores
- Pressure sore care



- | Stage 2: Scapula
- | Stage 3: Great trochanter
- | Stage 4: Sacrum bone
- | Unstageable: Heel

Perineal care

- Anatomically correct genitals for detailed care



Body positioning and transfer

- Communication
- Correct manual handling techniques
- Patient safety



- | Changing body positions
- | Supporting patient transfer

Medical Understanding

Airway Suction

- Understanding the procedure



- | Airway suction: oral, nasal and tracheotomy (*tube insertion only)
- | Tube feeding: oral, nasal and PEG (*tube insertion only)

Daily care with medical understandings



- | Stoma care
- | Bed bathing for the patients with central venous port
- | Perineal care for the patients with urinal catheters



Photo: Kyoto College of Nursing

* A wheelchair is not included.

DESCRIPTIONS

SPECIFICATIONS

Size: 149 cm / 58.7 inch Weight: 15kg / 33lbs

SET INCLUDES

1 manikin	1 male genitalia
1 denture	1 under pants
1 wig	1 wrench
1 oral cover for suctioning	1 instruction manual
1 female genitalia	

MATERIALS

Latex free
Soft resin

REPLACEMENT PARTS

11405-010	female genitalia
11405-020	male genitalia
11401-080	denture
11229-060	oral cover for suctioning
11229-050	lubricant for training models
11381-800	lubricant for catheterization

OPTIONAL PARTS

12746-020	tracheotomy tube
11251-320	simulated suppository set
11251-330	simulated feces set
11405-030	pajamas
11404-030	washable wig (elderly)
11229-070	sputum

RECOMMENDED DEVICES

Nasogastric tube feeding:
 -Nasal catheter 12fr -14fr
 -Gastrostomy catheter for replacements 20fr
Indwelling Urinary Catheter:
 16fr (male) / 14fr (female)