

MW2A



Male Catheterization/Enema Simulator

MW52



Female Catheterization/Enema Simulator II



MW2A



MW52

High quality simulators for urethral catheterization, perineal care and enema practice

Achieve better training results with our soft and realistic catheterization simulator that can also be used on SP or full body manikins.

FEATURES

- 1 | Soft and flexible genitalia allow proper cleaning and procedural preparation
- 2 | Successful insertion confirmation by fluid outflow
- 3 | Enema insertion in lateral position
- 4 | Detachable external genitalia


SKILLS

- | Intermittent catheterization
- | Indwelling catheterization
- | Perineal care
- | Enema
- | Localization of urethral meatus
- | Cleansing and aseptic technique
- | Catheter insertion and placement
- | Catheter Removal
- | Manual bladder compression

Compact sets consist of a genitalia unit and accessories for hybrid training and table top demonstration.

MW2C 11381-199-B

Male Genital Unit only



MW2D 11433-299-B

Female Genital Unit only



Enema Training



Pass of catheter can be observed through transparent bladder reservoir



DESCRIPTIONS

SPECIFICATIONS

Size: 34x33x24cm / 13.5x13x9.5 inch
Weight: MW2A 3.5kg / 7.7lbs,
MW52 3.4kg/7.5lbs

SET INCLUDES

- 1 lower torso manikin
- 1 genital organ unit
- 3 valves (one is incorporated in the assembled unit)
- 1 irrigation bottle

MATERIALS

Soft resin
Latex free

- 1 rectum tube
- 3 irrigation drainage tubes
- 1 lubricant
- 1 storage case
- 1 instruction manual

REPLACEMENT PARTS

- 11381-100-01 clip (small) for urethral tube for MW2A
- 11381-100-02 clip (large) for rectum unit for MW2A/2B/2C
- 11381-100-03 urinary bladder reservoir
- 11381-100-04 rectum unit
- 11381-300 6 replacement valves
- 11433-010 1 replacement female skin
- 11433-020 2 replacement female meatus
- 11381-800 lubricant
- 11381-410 1 replacement male skin
- 11381-510 2 replacement male meatus

RECOMMENDED DEVICES

* Catheters are not included.

