

MW32



NEW

# NCPR Simulator

Manikin that facilitates Cardio-Pulmonary Resuscitation training of newborn baby



## FEATURES

- 1 | Airway management training with observable bilateral/unilateral ventilation
- 2 | Soft and realistic umbilical cord capable of venous and arterial line placement, and blood sample collection
- 3 | Feeding tube placement confirmable through gastric fluid aspiration and
- 4 | Chest compression

## SKILLS

### Basic neonatal care

- | Changing body position
- | Hygiene
- | Palpation of the anterior fontanelle
- | Umbilical catheterization
- | Umbilical cord care
- | Tube feeding (nasal and oral) and gastric fluid aspiration
- | Clamping, tying and cutting of the umbilical cord
- | Suction (nasal, oral)

### Emergency Simulation

- | Chest compression
- | Endotracheal intubation
- | BVM ventilation

\*The product is not designed for bathing.

## ANATOMY

- | ribs
- | bronchi
- | lungs
- | stomach
- | esophagus
- | umbilical cord

## DESCRIPTIONS

### SPECIFICATIONS

Size: 50 cm / 19.6 inch Weight: 2.6kg / 5.7 lbs

### SET INCLUDES


- 1 newborn manikin
- 1 lubricant
- 1 instruction manual

### MATERIALS

Soft resin  
Latex free

### REPLACEMENT PARTS


- 11400-010 1 mask for MW21, MW32
- 11400-020 1 chest skin for MW21, MW32
- 11400-030 4 umbilical cords for MW21, MW32
- 11229-050 1 lubricant



Palpation of the anterior fontanelle



Umbilical cord care  
(one vein and two arteries)  
umbilical catheterization and blood sampling



Tube Feeding (nasal and oral)  
(tube placement confirmation with bubble sounds)



Airway suctioning (nasal and oral)



Chest compression



Endotracheal intubation  
(confirmation of unilateral/bilateral ventilation)



BVM ventilation  
(rise of the chest)

### RECOMMENDED DEVICES

Feeding tube: 5fr.  
Suction catheter: 5fr.  
Endotracheal tube: 3.0mm ID  
Umbilical catheter: 5fr., 30cm  
Laryngoscope: #0  
BVM: neonatal

\* Instruments are not included.

