


MW35



NEW

Cervical Assessment Simulator

Realistic yet durable material
for repetitive practice of cervix scoring



Station zero, the line that connects ischial spines can be recognized.

FEATURES

- 1 | Genital area made of innovative material that excels in elasticity, durability, and resilience, offers repetitive training.
- 2 | 5 types of cervical dilation models and fetal descent positions can be interchanged with a single action.
- 3 | Anterior/posterior fontanelles to determine the fetal position

SKILLS

- Cervix scoring in the first stages of labor
- Check fetal head position
- Determination of bishop score (cervix score)
- Speculum examination (i.e. Pap smear testing)
- Urinary catheterization
- Cervical cytology

KEY FEATURES



Cervical dilation insert can be attached by a single action.



Descent of the fetus can be set by the lever on the holder. (from -3 to +3)



External genitalia that excels in elasticity, durability and resilience

DESCRIPTIONS

SPECIFICATIONS

Size: W74xD44 x H33.5 cm
W29xD17xH13 inch

Weight: 10.3 kg/ 22.7 lbs

MATERIALS

Soft resin, hard plastic
Latex free

RECOMMENDED DEVICES

-

SET INCLUDES

1 maternal body	3 lubricant (150 ml)
1 genitalia unit	1 talcum powder
1 vagina unit	1 instruction manual
1 set of cervical dilation inserts (5 variations)	
1 holder for the cervical dilation inserts	

REPLACEMENT PARTS

11415-010	1 set of cervical dilation inserts (5 variations)
11415-020	3 lubricant (150 ml)
11415-030	1 genitalia unit (for pelvic examination and delivery assistance)

