

MW38



NEW

Uterine Involution Assessment Simulator

Practice highly realistic uterine involution assessments via a soft, durable and elastic postpartum torso



Palpation of uterine fundus



Measurement of fundal height

FEATURES


- 1 | Soft and elastic abdominal wall that permits highly realistic palpation to assess the degree of uterine involution and to measure fundal height.
- 2 | Four interchangeable uterine inserts

SKILLS

- | Assessment of uterus in early puerperium
- | Fundal height measurement
- | Palpation of uterine fundus
- | Perineal cleansing and assessment
- | Postpartum massage
- | Assessment of anal prolapse

KEY FEATURES

4 variations of uterine inserts




NORMAL (boundary between abdominal wall and uterine is hard and clear)

ABNORMAL (uterine is soft or unclear)

LANDMARKS

- | pubic symphysis
- | umbilicus
- | uterus



DESCRIPTIONS

SPECIFICATIONS

Size: W28 x D42 x H26 cm Weight: 7.4kg / 16.3 lbs
W11 x D16.5 x H10.2 inch

SET INCLUDES

- 1 | puerperal lower torso
- 1 | set of uterine inserts (4 variations)
- 1 | talcum powder
- 1 | instruction manual

MATERIALS

Soft resin
Latex free

RECOMMENDED DEVICES

-

REPLACEMENT PARTS

- 11418-010 1 skin cover (for uterine involution assessment)
- 11418-020 1 set of uterine inserts (4 variations)

