


MW3

Epidural Anesthesia Simulator

Practice and perform epidural anesthesia with confidence

Epidural Anesthesia Simulator offers multiple training methods on both thoracic and lumbar areas, including lumbar punctures.




FEATURES


- 1 | True-to-life resistance with the "pop" sensation
- 2 | Left and right lateral positions
- 3 | Simulates loss of resistance to saline or water
- 4 | The broad puncture pad allows paramedian access as well as median approach.
- 5 |

SKILLS


- | Epidural anesthesia procedures at lumbar and thoracic area
- | Lumbar puncture



Landmark palpation



Epidural anesthesia



Epidural anesthesia

DESCRIPTIONS

SPECIFICATIONS

Size: W48 x D23 x H29 cm
W18.9 x D9 x H11.4 inch

SET INCLUDES

- 1 | life size unisex torso
- 5 | puncture pads
- 1 | supporting stand
- 1 | 50ml syringe

Soft resin
Latex free

REPLACEMENT PARTS

11382-010 2 replacement puncture pads

RECOMMENDED DEVICES

Epidural needle: 18G
Spinal needle: 21G