

MW54/MW53/MW55


Strap-on Catheterization / Enema Simulator

Reality in skills training combined with scenarios for communication.



True-to-life softness

Labia can be kept open through the catheterization procedures with the non-dominant hand.



FEATURES

- 1 | Soft female genitalia permits holding the labia in open position throughout the procedure from sterilization to the placement of the catheter.
- 2 | Anatomically correct structure facilitates acquisition of body sense along with the procedure.
- 3 | Compatible with balloon catheters. Water (simulated urine) flows out when appropriate length of a catheter is inserted.

SKILLS

- | Cleansing and aseptic technique
- | Catheter insertion and placement
- | Confirmation of urine discharge
- | Catheter removal
- | Enema



The catheter can be inserted when the penis is held at the right angle.



Injection of enema fluid



Realistic resistance to the catheter

DESCRIPTIONS

SPECIFICATIONS

Life-sized genitalia model: W8 x D21 x H18cm, 0.5kg (male)
W8 x D21 x H9cm, 0.4kg (female)

SET VARIATIONS

female: 11434-000
male: 11435-000
value set: 11436-000

REPLACEMENT PARTS

11381-800 lubricant

SET INCLUDES

1 genitalia model	1 lubricant
1 reservoir bag	1 catheter for maintenance
1 drainage tube	1 instruction manual
1 strap band	
1 shorts	