

MW58

Male Suprapubic Catheter Care Simulator

NEW

MW59

Female Suprapubic Catheter Care Simulator

NEW

A comprehensive trainer for suprapubic catheter changing.



FEATURES

- 1 | Changing suprapubic catheters
Successful insertion confirmation by fluid outflow
Training includes check points in the skills such as catheter mobility and confirmation of urine discharge.
- 2 | Suprapubic catheter care (securing catheter in place)
Facilitates training of securing catheter with gauze and surgical tape.
Confirmation of catheter placement
- 3 | Observation of the catheter through transparent bladder wall.

SKILLS

- | Changing suprapubic catheters
- | Suprapubic catheter care (securing catheter in place)
- | Placement of transurethral balloon catheters
- | Enema at lateral positions

MW59




MW58



Confirmation of discharge of urine.



Catheter insertion



DESCRIPTIONS

SET INCLUDES

- 1 lower torso manikin
- 1 genital organ unit
- 4 valves (2 are incorporated in the assembled unit)
- 1 irrigation bottle
- 3 irrigation drainage tubes
- 1 instruction manual

REPLACEMENT PARTS

- 11440-010 SPC male skin
- 11441-010 SPC female skin
- 11441-020 cap for suprapubic opening
- 11441-030 SPC urinary bladder reservoir