

M67C

11297-300

## Strap-on Intramuscular Injection Trainer

For practice between students, simulated patients and mannikin, facilitating their understanding of the intramuscular injection process and communication



### FEATURES

- 1 | The water tank is located inside this trainer
- 2 | Realistic tactile feeling of skin and acromion allows palpating the anatomical landmark
- 3 | Easy to drain

### SKILLS

- | Intramuscular Injection training
- | COVID-19 vaccination training

### DESCRIPTIONS

#### SPECIFICATIONS

Size: W10.5 x D26 x H5.5 cm  
W4.2 x D10.3 x H2.2 in.

#### SET INCLUDES

- 1 main unit  
(frame, skin, soft tissue, muscle,  
acromion, drainage reservoir, belts)
- 1 instruction manual

#### MATERIALS

Soft resin  
Latex free

#### RECOMMENDED DEVICES

IM needle 23G

#### REPLACEMENT PARTS

- 11297-220 replacement skin
- 11297-230 replacement soft tissue part
- 11297-240 replacement muscular part

MW64

11454-989,999

## Vaccination Training Set

MW64 Vaccination Training Set will undeniably help in training thousands of healthcare workers who are involved in COVID-19 vaccination. It's simple, easy to use and low-maintenance



### FEATURES

- 1 | The landmark acromion process is palpable, allowing for realistic practice of IM injection
- 2 | The injection pad has a highly absorbent sponge built in, so that you can actually inject the liquid repeatedly
- 3 | This vaccination training set can be used together with Kyoto Kagaku patient care manikins, SPs and trainees

### SKILLS

- | Deltoid intramuscular injection
- \*The pad's depth is 1 1/2" (38.1mm)

### DESCRIPTIONS

#### SPECIFICATIONS

Size: W12 x D14 x H5.5 cm  
W4.7 x D5.5 x H2.2 in.

#### SET INCLUDES

- 5 vaccination training pads
- 5 pad holders
- 5 strap-on belts
- 1 instruction manual

#### MATERIALS

Soft resin  
Latex free

#### RECOMMENDED DEVICES

IM needle 23G

#### REPLACEMENT PARTS

- 11454-010 Pad for Vaccination Training Set (Set of 5)
- 11454-020 Pad for Vaccination Training Set